



Clinical Collaboration Platform

**Vue PACS and Vue Archive 12.2
IHE Integration Statement**

**Part # AD0551
2018-04-10**

Integrating the Healthcare Enterprise (IHE) Integration Statement			
Vendor	Product Name	Version	Date
Carestream Health, Inc	CARESTREAM Vue Archive (Vue for Clinical Collaboration Platform), CARESTREAM Vue PACS: Client Suite, Workflow Manager, IMS, CD Direct Suite, XDS Registry.	CARESTREAM Vue PACS, Vue Archive v12.2	2018-04-10
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed in this table.			
Note: This document is intended for use by Carestream customers.			

Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow (SWF)	Image Manager/Image Archive Image Display PPS Manager	None
Imaging Object Change Management (IOCM)	Image Manager/Image Archive Change Requestor	None
Access to Radiology Information (ARI)	Image Manager/Image Archive Image Display Report Repository External Report Repository Access	None
Patient Information Reconciliation (PIR)	Image Manager/Image Archive PPS Manager	None
Consistent Presentation of Images (CPI)	Image Manager/Image Archive Image Display	None
Cross-enterprise Document Sharing (XDS.b)	Document Consumer Document Repository Document Registry Document Source	Multiple document
Cross-enterprise Document Sharing for Imaging (XDS.b-l)	Imaging Document Consumer Imaging Document Source	Set of DICOM Instances, Multiple documents None
Cross-enterprise Document Reliable Interchange (XDR)	Document Recipient Document Source	Accept Limited Metadata
Cross-Enterprise Sharing of Scanned Documents (XDS-SD)	Content Consumer	None
Basic Patient Privacy Consents (BPPC)	Content Consumer	None

Integration Profiles Implemented	Actors Implemented	Options Implemented
Patient Demographic Query HL7 V3 (PDQv3)	Patient Demographics Consumer	None
Patient Identifier Cross-referencing (PIX)	Patient Identifier Cross-reference Consumer	None
Evidence Documents (ED)	Image Display Image Manager/Image Archive	None
Cardiology Evidence Documents (ED-CARD)	Image Display Image Manager/Image Archive	None
Simple Image and Numeric Report (SINR)	Report Repository External Report Repository Access	None
Mammography Image (MAMMO)	Image Manager/Image Archive Image Display	None
Key Image Note (KIN)	Evidence Creator Image Display Image Manager/Image Archive	None
Nuclear Medicine Imaging (NMI)	Image Manager/Image Archive	None
Audit Trail and Node Authentication (ATNA)	Secured Application Secured Node	None
Consistent Time (CT)	Time Client	None
Portable Data for Imaging (PDI)	Portable Media Creator Portable Media Importer Image Display Display	Web Content None None None
Reporting Workflow Profile (RWF)	Image Manager/Image Archive	None
Import Reconciliation Workflow Integration Profile (IRWF)	Image Manager/Image Archive	None
Teaching File and Clinical Trial Export (TCE)	Receiver	None
Post-Processing Workflow (PWF)	Image Manager/Image Archive	None
Presentation of Grouped Procedures Integration Profile (PGP)	Image Manager/Image Archive	None

Link to Carestream Health IHE information:		http://ihe.carestream.com
Links to Standards Conformance Statements for the Implementation		
DICOM	http://dicom.carestream.com	
Links to general information on IHE		

Link to Carestream Health IHE information: http://ihe.carestream.com		
Links to Standards Conformance Statements for the Implementation		
In North America: https://www.rsna.org/ihe.aspx	In Europe: https://www.ihe.net	In Japan: www.jira-net.or.jp/ihe-j



Authorized Representative (European Union)



Carestream Health France
1, rue Galilée
93192 NOISY-LE-GRAND CEDEX
FRANCE



Importer for European Union

Carestream Health Netherlands B.V.
Bramenberg 12
3755 BZ Eemnes
The Netherlands



150 Verona Street
Rochester, New York, USA 14608

© Carestream Health, Inc., 2018

CARESTREAM is a trademark of Carestream Health

Made in the USA